# Quick Setup Guide Transmitter M200

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---

### Content

1	Safety and Intended Use	31
2	Technical Data and Further Information	32
3	Mounting	32
4	Electrical Connection	32
5	Menu Structure	34
6	Operation On-Site	35
7	Sensor Calibration	37
8	Maintenance	38
9	Disposal	38
10	EC-Declaration of Conformity	38

# 1 Safety and Intended Use

 **Note:** The Quick Setup Guide is a brief operating instruction.

The M200 transmitter must only be installed, connected, commissioned, and maintained by qualified specialists e.g. electrical technicians in full compliance with the instructions in this Quick Setup Guide, the applicable norms and legal regulations.

The specialist must have read and understood this Quick Setup Guide and must follow the instructions it contains. If you are unclear on anything in this Quick Setup Guide, you must read the Operation Manual (supplied on CD-ROM). The Operation Manual provides detailed information on the instrument.

The M200 transmitter should be operated only by personnel familiar with the transmitter and who are qualified for such work.

## Intended Use

The M200 is a 4-wire transmitter for analytical measurements with 4 (0) to 20 mA output signal. The M200 is a multi-parameter transmitter for pH/ORP, conductivity, dissolved oxygen and ozone measurement. It is available as 1-channel or 2-channel version and compatible with easySense™ and ISM® sensors.

The M200 transmitter is designed for use in the process industries, in non-hazardous areas.

## M200 parameter fit guide

Parameter	ISM	easySense
pH/ORP	•	pH 31, pH 32, pH 33, pH 34, ORP 41
pH/pNa	•	–
Conductivity 2-e	•	Cond 71, Cond 72, Cond 73
Conductivity 4-e	•	Cond 77
Amp. Dissolved Oxygen ppm	•	O <sub>2</sub> 21
Dissolved Ozone	•	–

## 2 Technical Data and Further Information

The most important technical data such as supply voltage are given on the name plate either outside or inside the transmitter housing. For further technical data such as accuracy refer to the Operation Manual. This document, the Operation Manual and software are on the supplied CD-ROM. You can also download the documentation via Internet "www.mt.com/M200".

## 3 Mounting

The M200 transmitter is available as ½ DIN and ¼ DIN version. For installation drawings refer to the Operation Manual.

1. ½ DIN only: Mount the supplied cable glands at the housing.
2. Mount the transmitter. You have the following possibilities:
  - Panel mounting: ½ DIN and ¼ DIN version
  - Wall mounting: ½ DIN
  - Pipe mounting: ½ DIN

## 4 Electrical Connection

**DANGER! Mortal danger by electric shock:** Power off instrument during electrical connection.


**Note:** This is a 4-wire-product with an active 4–20 mA analog output. Do not supply power to terminal no. 1 to 6 of terminal block TB2.

1. Switch off supply voltage.
2. Connect mains supply as follows:
  - 20 to 30 V DC: **-N** for Neutral and **+L** for Line
  - 100 to 240 V AC **-N** for Neutral and **+L** for Line
3. Connect relay output signals, analog output signals and digital input signals according to the Operation Manual.
4. Connect sensor signal to terminal block **TB3** and **TB4** (only 2-channel version) according to the Operation Manual.


## Terminal block (TB) definitions

### ½ DIN housing (wall mount)


- 1 USB-connection for software update
- 2 TB1: Relay outputs
- 3 Power connection
- 4 TB2: Analog outputs and digital inputs
- 5 TB3: Sensor connection
- 6 TB4: Sensor connection (only 2-channel version)

### ¼ DIN housing (panel mount)


- 1 Power connection
- 2 TB4: Sensor connection (only 2-channel version)
- 3 TB3: Sensor connection
- 4 TB2: Analog outputs and digital inputs
- 5 TB1A, TB1B: Relay outputs
- 6 USB-connection for software update

## 5 Menu Structure


## 6 Operation On-Site

Left: Measurement mode (example), Right: Edit mode (example)

- 1 Channel information  
A: 1st connected ISM or easySense sensor.  
B: 2nd connected ISM or easySense sensor (only 2-channel version).  
H: Transmitter is in Hold mode. Refer to the Operation Manual of the transmitter.
- 2 1st line (a), standard configuration
- 3 2nd line (b), standard configuration
- 4 3rd line (c), Measurement mode: The screen depends on configuration.  
Edit mode: Navigation through the menu or editing parameters
- 5 4th line (d): Measurement mode: The screen depends on configuration.  
Edit mode: Navigation through the menu or editing parameters
- 6 If a  $\uparrow$  is displayed, you can use the  $\blacktriangleright$  or the  $\blacktriangleleft$  key to navigate to it. If you click [ENTER] you will navigate backwards through the menu (go back one screen).

You can configure the information shown on the display for each line. By default in the measurement mode no values are shown in the third and in the fourth line of the display. For configuration refer to the Operation Manual of the transmitter.

Navigation keys	Description
	<ul style="list-style-type: none"> <li>• Enter Menu mode.</li> <li>• Navigate backwards within a changeable data entry field.</li> </ul>
	<ul style="list-style-type: none"> <li>• Enter Calibration mode.</li> <li>• Navigate forward within a changeable data entry field.</li> </ul>
	<p>Return to Measurement mode. Press the ◀ and ▶ key simultaneously (escape).</p> <p><b>NOTE:</b> In order to back up only one menu page, move the cursor under the UP Arrow character (↑) at the bottom right of the display screen and press [ENTER].</p>
	<ul style="list-style-type: none"> <li>• Increase a digit.</li> <li>• Navigate within a selection of values or options of a data entry field.</li> </ul>
	<ul style="list-style-type: none"> <li>• Enter Info mode.</li> <li>• Decrease a digit.</li> <li>• Navigate within a selection of values or options of a data entry field.</li> </ul>
	Confirm action or selection.

**Note:** Some screens require configuring multiple values via the same data field (ex: configuring multiple set points). Be sure to use the ▶ or ◀ key to return to the primary field and the ▲ or ▼ key to toggle between all configuration options before entering to the next display screen.

"Save changes" dialog	Description
Yes & Exit	Save changes and exit to measurement mode
Yes & ↑	Save changes and go back one screen
No & Exit	Do not save changes and exit to measurement mode

## 7 Sensor Calibration

 **Note:** For best process calibration results observe the following points. Take grab sample as close as possible to the measurement point of the sensor. Measure the sample at process temperature.

For a detailed description of the “Process”, “1-Point” and “2-Point” calibration method see Operation Manual of the M200 transmitter. Via iSense software you can calibrate the sensor with the “1-Point” or “2-Point” calibrate method. See Operation Manual of the iSense software.

As soon as the calibration is in progress no other calibration can be started.

### **Sensor Calibration Menu**

After every successful calibration the options “Yes” or “No” are available. After selection one of these options the message “Re-install sensor and Press [ENTER]” appears.

After pressing [ENTER] the M200 returns to the Measurement mode.

- **Yes**

Calibration values are stored in the sensor and used for the measurement.

- **No**

Calibration values are discarded.

## **8 Maintenance**

The transmitter requires no maintenance.

Clean the surfaces with a soft damp cloth and dry the surfaces with a cloth carefully.

## **9 Disposal**

Observe the applicable local or national regulations concerning the disposal of "Waste electrical and electronic equipment".

Disassemble the transmitter according to resources. Sort resources and supply them to recycling. Non-recyclable materials have to be disposed of in an environmental-friendly manner.

## **10 EC-Declaration of Conformity**

The EC-Declaration of conformity is part of the delivery.

ISM is a registered trademark of the Mettler-Toledo Group in Switzerland, Brazil, USA, China, European Union, South Korea, Russia and Singapore.

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© Mettler-Toledo GmbH 10/2019.
**Mettler-Toledo GmbH**
**Mettler-Toledo Thornton, Inc.**