METTLER TOLEDO

3201-UPW & 3201i-UPW pH Combination Electrode

Instruction Manual

cs	Návod k obsluze	2
da	Brugsanvisning	6
de	Betriebsanleitung	10
en	Instruction manual	14
es	Manual de instrucciones	18
fi	Käyttöohje	22
fr	Instructions d'utilisation	26
hu	Használati utasítás	30
it	Istruzioni per l'uso	34
jp	取扱説明書	38
ko	사용 매뉴얼	42
nl	Gebruikershandleiding	46
pl	Instrukcja obsługi	50
pt	Manual de instruções	54
ru	Инструкция по эксплуатации	58
sv	Bruksanvisning	62
th	คู่มือคำแนะนำ	66
zh	说明书	70

3201-UPW a 3201 i

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
Instruction Manual

Contents

1	Introduction	15
2	Safety instructions	15
3	Product description.....	15
4	Installation and commissioning	16
5	Operation.....	16
5.1	Calibration of the sensor and pH transmitter.....	16
5.2	Temperature sensor	16
6	Maintenance	16
7	Trouble-shooting	17
8	Environmental protection.....	17
9	Warranty	17

ISM is a registered trademark of the METTLER TOLEDO Group in Switzerland, USA, the European Union and a further five countries.
 pHure Sensor is a registered trademark of the METTLER TOLEDO Group in the USA and China.

1 Introduction

Mettler-Toledo Thornton pH sensors type 3201-UPW are low-maintenance electrodes with pressurized gel-electrolyte, silver-ion trap and built-in temperature sensor.

2 Safety instructions

The sensor consists primarily of glass and should be handled carefully. If the sensor is to be cleaned or calibrated using acid or alkaline solutions, safety goggles and gloves should be worn. Local safety regulations are applicable to the handling of the sensor. Please take extra care during installation/fitting of the sensor to ensure that the polymer thread (Pg 13.5) does not suffer any damage.

We recommend that the sensor may only be in operation with the original accessories from Mettler-Toledo Thornton. Operation and maintenance should be carried out only by trained personnel and staff familiar with the unit and who have read and understood this instruction manual.

3 Product description

The markings on each sensor provide the following information:

METTLER TOLEDO THORNTON	Sensor manufacturer
3201-UPW/i/120/Pt 1000	Type designation / i = ISM digital/shaft length in mm/temperature sensor
combination pH	Type of sensor
pH 1...11	pH measuring range
0...80 °C (176 °F)	Temperature range for operation
Order No. 52 00X XXX	Article number for ordering

Additionally, the sensor is supplied with a serial number on the red connector part for identification of each sensor produced.

4 Installation and commissioning

1. On unpacking, check the sensor for possible mechanical damage of the pH sensitive glass, the ceramic reference junction/diaphragm, the shaft and the VarioPin (VP) or ISM connector. Report any damage immediately to your Mettler-Toledo Thornton supplier.
2. Remove the watering cap and briefly rinse the sensor with deionized water. Important: Carefully and thoroughly remove the protective droplet of silicone sealant on the reference junction/diaphragm with the enclosed knife, as shown on the green pictures on the packing.
3. Again rinse the sensor well with deionized water. After rinsing, the sensor should only be dabbed dry with a tissue. Do not rub the pH sensitive glass, since this can lead to electrostatic charging and sluggish response times.
4. Carefully check the space behind the pH sensitive glass for the presence of any air bubbles and remove same by gently swinging the sensor in a vertical plane.
5. Please refer to the instruction manual for the relevant housing in order to install the sensor correctly.
6. Connect the sensor to the pH preamp or transmitter, using a cable with VarioPin (VP) or ISM connector. Please refer to the wiring and terminal diagram accompanying the cable.

5 Operation

5.1 Calibration of the sensor and pH transmitter

For 3201-UPW sensors, a 2-point calibration is recommended. After removing the watering cap, the sensor is dipped in succession into two different buffer solutions with known pH reference values and the transmitter is calibrated to these buffer values.

5.2 Temperature sensor

The integrated Pt1000 temperature sensor is primarily intended for automatic temperature compensation of the pH signal during calibration and normal operation. For accurate indication or control of the process temperature, a further sensor temperature calibration may be performed.

6 Maintenance

1. Carefully rinse the sensor tip and reference junction/diaphragm with deionized water after removal from the process. Under no circumstances should measuring solution containing solids be allowed to dry on these parts!
2. When the sensor is not in operation, store it with sensor tip and reference junction/diaphragm well submerged into 3M KCl.
3. If the sensor is stored mounted in its housing, the same rule as in 2) above applies, but the storage electrolyte has to be slightly modified, with 2 parts

of buffer solution pH 9 to 10 parts of electrolyte, this in order to avoid corrosion of housing parts.

4. If a sensor is stored dry for a few days by mistake, let it soak in the normal storage electrolyte for several hours before use.
5. Occasionally check the connector for possible traces of moisture. If necessary, clean well with deionized water or alcohol, and afterwards dry carefully.

7 Trouble-shooting

Sluggish response time (e.g. after a long and dry storage period) is an indication that reactivation of the pH sensitive glass may be necessary. Dip the pH sensitive glass (but no further sensor parts!) into Reconditioning (HF) Solution (METTLER TOLEDO product, article no. 51 319 001) for 5 min. Afterwards, rinse well with deionized water and let sensor tip and diaphragm soak in reference electrolyte for minimum 12 hours.

8 Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

9 Warranty

In the event of manufacturing faults, 12 months warranty from date of production is granted.

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05/2018 © METTLER TOLEDO. All rights reserved
- France** **Mettler-Toledo Analyse Industrielle S.A.S.**
- Germany** **Mettler-Toledo GmbH**
- Switzerland** **Mettler-Toledo (Schweiz) GmbH**
Mettler-Toledo GmbH, Process Analytics, Im Hackacker 15, CH - 8902 Urdorf