


# Onsite Service

## For All Pipette Brands...

**Rainin, Gilson, Eppendorf, Thermo Fisher, Sartorius...**  
**Your pipettes need regular check ups. Whatever your service need, trust Rainin for exceptional quality and attention to detail.**


Rainin makes it easy to get fast, reliable pipette service and calibration.

- We come to you
- Eliminate chain-of-custody concerns
- Nationwide service capability
- ISO 17025 accredited by A2LA
- Calibrated to ISO 8655 or Rainin on-site standard tolerances
- METTLER TOLEDO balances
- Minimal downtime
- Calibration in the environment where pipettes are used



Service Plan

Calibration Plans



See reverse page for detail information

# RAININ

Feel the difference!

## Ordering Information

BASIC PLAN	PREVENTIVE MAINTENANCE	AS FOUND	AS RETURNED	CALIBRATION CERTIFICATES
When measurement data is not required	✓	n/a	4 test weighings @ 10% and 100%	N/A

- Full preventive maintenance (PM) with seal replacement and leak test
- Calibrated to Rainin on-site standard tolerances

ORDER No.	DESCRIPTION	PRICE
<b>30445525</b>	Single-channel	\$ 28.80
<b>30445526</b>	Multichannel All – 8 channels	55.40
<b>30445527</b>	Multichannel All – 12 channels	68.60

STANDARD PLAN	PREVENTIVE MAINTENANCE	AS FOUND	AS RETURNED	CALIBRATION CERTIFICATES
For labs with routine service needs	✓	n/a	4 test weighings @ 10% and 100%	Standard Certificate

- Full preventive maintenance (PM) with seal replacement and leak test
- Calibrated to ISO 8655 or Rainin on-site standard tolerances
- Calibration certificate with detailed measurement data

ORDER No.	DESCRIPTION	PRICE
<b>17800347</b>	Single-channel	\$ 39.00
<b>17800292</b>	Multichannel All – 8 channels	108.40
<b>17800256</b>	Multichannel All – 12 channels	126.00

ADVANCED PLAN	PREVENTIVE MAINTENANCE	AS FOUND	AS RETURNED	CALIBRATION CERTIFICATES
When prior performance matters	✓	4 test weighings @ 10% and 100%	4 test weighings @ 10% and 100%	Standard Certificate

This plan includes all services provided under the Standard Plan.

Plus:

- As found calibration results – including failure confirmation

ORDER No.	DESCRIPTION	PRICE
<b>17800320</b>	Single-channel	\$ 47.20
<b>17800272</b>	Multichannel All – 8 channels	145.60
<b>17800236</b>	Multichannel All – 12 channels	177.60

ACCREDITED PLAN	PREVENTIVE MAINTENANCE	AS FOUND	AS RETURNED	CALIBRATION CERTIFICATES
For complete ISO 17025 documentation	✓	4 test weighings @ 10%, 50% and 100%	4 test weighings @ 10%, 50% and 100%	ISO 17025 Certificate

This plan includes all services provided under the Advanced Plan.

Plus:

- Calibrated to ISO 8655 tolerances
- ISO 17025 calibration certificate incl. uncertainty of measurement (certified by A2LA), historical trending and Calibration label

ORDER No.	DESCRIPTION	PRICE
<b>30451900</b>	Single-channel	\$ 86.60
<b>30451876</b>	Multichannel All – 8 channels	241.60
<b>30451852</b>	Multichannel All – 12 channels	302.40

**Prices do not include applicable travel fees.**

For more information on service selection, please visit [www.RaininService.com](http://www.RaininService.com)

Phone orders: **1 800 4 RAININ (1 800 472 4646)**